

Luxator (Thin Blade)

By cutting the periodontal ligament, easily dislocate from the root. Advantage in saving form of a skeleton and minimizing the damage on the tissue.



MLB-03 Luxating Elevator B1 (Right)

- Blade width 3,0mm
- For posterior



MLB-01 Luxating Elevator B1 (Straight)

- Blade width 3,0mm
- For anterior



MLB-02 Luxating Elevator B1 (Left)

- Blade width 3,0mm
- For posterior



MLB-04 Luxating Elevator B1 (Curved)

- Blade width 3,0mm
- For edge in posterior area