

Luxator (Bold Blade)

By cutting the periodontal ligament, easily dislocate from the root. Advantage in saving form of a skeleton and minimizing the damage on the tissue.



MLC-03 Luxating Elevator C1 (Right)

- Blade width 5.0mm
- For posterior



MLC-01 Luxating Elevator C1 (Straight)

- Blade width 5.0mm
- For anterior



MLC-02 Luxating Elevator C1 (Left)

- Blade width 5.0mm
- For posterior



MLC-04 Luxating Elevator C1 (Curved)

- Blade width 5.0mm
- For edge in posterior area