

# Elevator

Use to provoke in start of extraction, or to lift and remove from bone



**MEA-06** Elevator A6

- Width 4.9mm
- For Upper and lower Posterior and roots



**MEA-07** Elevator A7

- Width 2.7mm
- For deep tooth roots



Curved, Backward



**MEA-08** Elevator A8

- Width 3.5mm
- Apex and small tooth in Upper Posterior area



Curved, Backward



**MEA-09** Elevator A9

- Width 4.3mm
- For Upper and lower Canine



**MEA-10** Elevator A10

- Width 4.3mm
- For Upper and lower Canine

