

# Elevator

Use to provoke in start of extraction, or to lift and remove from bone



**MEB-01** Elevator B1

- Width 3,0mm
- Use for fractured tooth during extraction or to remove deep tooth roots



**MEB-02** Elevator B2

- Width 3,0mm
- Use for fractured tooth during extraction or to remove deep tooth roots



**MEC-01** Elevator C1

- Width 3,0mm
- Use for fractured tooth during extraction or to remove deep tooth roots



**MEC-02** Elevator C2

- Width 3,0mm
- Use for fractured tooth during extraction or to remove deep tooth roots